

Keys	Function	Comments
K	K-complex ruler	
L	75µV bars	
W	score "Wake"	
1	score "N1"	
2	score "N2"	
3	score "N3"	
R	score "REM"	
M	score "Movement Arousal"	
A	mark epoch as artefact	
C	reset current epoch	
H	show hypnogram	
E	show polysomnography statistics	
up	zoom all +	
down	zoom all -	
left	go back one epoch	
right	go forward one epoch	
page up	go forward 10 epochs (fast forward)	
page down	go back 10 epochs (fast backward)	
Pos 1	go to epoch 1	
end	go to last epoch	
ctrl+s	save as/save	
ctrl+i	import hypnogram (.txt)	
ctrl+e	export hypnogram (.txt)	
F6	show cluster plot	
F7	display filters on/off	
F8	dc remove on/off	
F9	highlight K-complexes on/off	
F12	show electrode configurator	select/deselect individual electrodes, set zoom for each electrode